



Picture 1 – Unlicensed Consumable. – waveform fails to achieve Biphasic breath for set iTIME.

Picture 2 Licensed Carefusion Consumable. Biphasic breath delivered to patient for set iTIME.



Picture 1 – unlicensed consumable – difference of 1.4cm is a failed Biphasic breath.

Picture 2 – licensed consumable – difference of 3.1cm is a biphasic breath reaching the patients alveolar.

The evidence above shows that the use of unlicensed consumables limits the SiPAPs ability to deliver optimal ventilation to your patient.

The use of unlicensed consumables has been shown to deliver unstable pressures to the patient and potentially lengthen patients stay on the unit.

A live demonstration of the above has been offered to me which I have:

Accepted

Declined

Having being provided with the information above it is our decision to continue to use the unlicensed consumables.

Signed on behalf of the trust: _____ Print name _____

Date _____ Position in the trust _____